


# Reaction Optimization Using a Microwave Autosampler

## Introduction

After identifying a novel reaction, the time consuming and often tedious process of reaction optimization begins, followed by substrate screening. While microscale reactions can be performed in parallel using a heating block, some reaction parameters cannot be adjusted easily in this format, consuming precious time for manually performed, sequential reactions. To overcome this process the CEM Explorer auto-sampler was developed, allowing researchers a more efficient way to survey and improve their chemistry.

## Procedure

The Hantzsch dihydropyridine synthesis (Scheme 1) is a common multicomponent reaction in which ammonia, an aldehyde, and (most often) a  $\beta$ -keto ester undergo a series of condensations in a single flask. The resultant 1,4-dihydropyridine compound is frequently isolated, however, it can spontaneously oxidize to the corresponding substituted pyridine.



**Scheme 1. General scheme for the Hantzsch dihydropyridine synthesis**

Beginning with 28% ammonium hydroxide solution (1.0 equiv), ethyl acetoacetate (2.0 equiv), and benzaldehyde (1.0 equiv), the reaction was performed according to literature, general microwave conditions<sup>1</sup>. In our hands, the method produced erratic results (Table 1; entries 1 and 2) of inconsistent temperatures and poor conversions to product (significant starting material remained and only minor product was formed in both entries as visualized by TLC [30% EtOAc/Hex]). Since the initial reactions resulted in an elevated pressure (214-260 psi), lower temperatures were surveyed (entries 3-5). A decrease in the isolated yield of desired product was noted with higher temperature, corresponding to an increase in unidentified byproducts. Performing the reaction at 150 °C resulted in an 83% isolated yield. At lower temperatures, a significant amount of starting materials remained while the product yield dropped dramatically. A longer reaction time at 150 °C did not result in additional product formation, while incomplete conversion was seen with shorter reaction times (entries 6 and 7, respectively). The use of ethanol as a solvent (entry 8) or excess ammonia (entry 9) did not improve the amount of isolated product, so previous conditions were used (150 °C, 5 minutes). After a rapid optimization, the product was furnished in an 83% isolated yield (entry 4). While this series of reactions was performed without automation, it is easy to adapt a reaction plan of nine entries to the Explorer platform, saving time and energy.


**Table 1. Optimization of microwave assisted Hantzsch dihydropyridine synthesis**

Entry	Temperature [°C]	Time [mm:ss]	Yield [%] <sup>a</sup>
1	198 <sup>[b]</sup>	1:40	—
2	188 <sup>[b]</sup>	1:40	—
3	170	5 <sup>[c]</sup>	68
4	150	5 <sup>[c]</sup>	83
5	130	5 <sup>[c]</sup>	53
6	150	10 <sup>[c]</sup>	83
7	150	3 <sup>[c]</sup>	70
8	150	5 <sup>[c]</sup>	73
9	150	5 <sup>[c]</sup>	59


<sup>[a]</sup> Isolated yield crude product <sup>[b]</sup> Reaction conditions: Fixed Power 45 W; 250 °C safe temperature; 300 psi; total time 100 s

<sup>[c]</sup> Hold time

Once optimal conditions were determined (Standard Method; 150 °C; 5 min hold time), varied electronic and steric substrates were selected. To perform a small library synthesis, four 10-mL microwave reaction vials were prepared by adding the correct amounts of the specified reagents (Table 2; entries 2-5), and a stir bar to each vessel. All vials were capped and placed in an Explorer rack, then programmed with the optimized method.



**Table 2. Substrate scope of neat, microwave assisted Hantzsch dihydropyridine synthesis**

Entry	R <sup>1</sup>	R <sup>2</sup>	R <sup>3</sup>	Yield [%] <sup>a</sup>
1	OEt	Me	Ph	83
2			Ph	99
3	OEt	Me	<i>p</i> -MeOC <sub>6</sub> H <sub>4</sub>	68 <sup>[b]</sup>
4	OEt	Me	2-furyl	86
5	OEt	Me	2-pyridyl	60

<sup>[a]</sup> Isolated yield crude product <sup>[b]</sup> Reaction required 10 minutes

Varied 1,4-dihydropyridines were synthesized in moderate to good yields in a “hands-free” fashion, using the Explorer auto-sampler to run this small library of compounds in under an hour. Using the CEM Explorer 96, up to ninety-six combinations of starting materials can be prepared and screened in one session, making this instrument ideal for methodology, combinatorial chemistry, and reaction optimization.

<sup>†</sup>Torchy, S.; Cordonnier, G.; Barbry, D.; Vanden Eynde, J. J. Hydrogen Transfer from Hantzsch 1,4-Dihydropyridines to Carbon-Carbon Double Bonds under Microwave Irradiation. *Molecules*. 2002, 7, 528–533.



**United States (Headquarters):**

CEM Corporation  
PO Box 200  
Matthews, NC 28106  
800-726-3331  
704-821-7015  
Fax: 704-821-7894  
info@cem.com  
cem.com

**France:**

CEM µWave S.A.S.  
Immeuble Ariane  
Domaine Technologique de Saclay  
4, rue Rene' Razel, 91892 ORSAY Cedex  
33 (01) 69 35 57 80  
Fax: 33 (01) 60 19 64 91  
info.fr@cem.com

**Germany, Austria, Switzerland:**

CEM GmbH  
Carl-Friedrich-Gauss-Str.9,  
47475 Kamp-Lintfort  
(49) 2842-9644-0  
Fax: (49) 2842-9644-11  
info@cem.de

**Ireland:**

CEM Technology (Ireland) Ltd.  
Sky Business Centre,  
9a Plato Business Park,  
Damastown, Dublin 15  
+353 (0) 1 885 1752  
Fax: +353 (0) 1 885 1601  
info.ireland@cem.com

**Italy:**

CEM S.R.L.  
Via Dell' Artigianato, 6/8  
24055 Cologno al Serio (Bg)  
(39) 35-896224  
Fax: (39) 35-891661  
info.srl@cem.com

**Japan:**

CEM Japan K.K.  
2-18-10 Takanawa, Minato-ku, Tokyo  
108-0074  
+81-3-5793-8542  
Fax: +81-3-5793-8543  
info@cemjapan.co.jp

**United Kingdom:**

CEM Microwave Technology Ltd.  
2 Middle Slade, Buckingham Industrial  
Estate, Buckingham MK181WA  
(44) 1280-822873  
Fax: (44) 1280-822873  
info.uk@cem.com