



## Cl-MPA ProTide and Cl-TCP(Cl) ProTide Resin Loading Procedure

This procedure details the simple and efficient method of loading Cl-MPA ProTide and Cl-TCP(Cl) ProTide resins with the Liberty series automated microwave peptide synthesizers at the 0.10 mmol scale. For the Parameter Values for other synthesis scales, contact [synthesis.support@cem.com](mailto:synthesis.support@cem.com).

This document should be used in conjunction with the Liberty Manual and the Safety Data Sheet (P/N 601400). Read and fully understand all documentation before operating the instrument.




Cl-MPA ProTide



Cl-TCP(Cl) ProTide

### ⚠ WARNING

Proper precautions must be taken to avoid contact with reagents or reagent vapors. Protective gear should be worn as outlined in the user's safety program for hazardous materials and the reagent manufacturer's safety data sheet. Refer to these guidelines for proper handling and disposal of the reagents. Dispose of all waste in accordance with all applicable local, state, and federal health and safety




## Step One: Program Cycle

Create or edit the “0.10 Chloride Loading (HS)” cycle editor.

1. From the “Edit” tab, select “Cycles...”
2. Select the following filters: “0.10” synthesis scale, “Amino Acid” and “High-Swelling”.
3. Highlight the “Chloride Loading (HS)” or “Trityl DCA Loading (HS)” cycle. If applicable, rename “Trityl DCA Loading (HS)” to “Chloride Loading (HS)”. NOTE: If using software v1.5 or earlier, the Coupling Cycle Step will read “Coupling (Trityl Loading)” and the Microwave Method Parameter Value will read “DCA Trityl Loading”.
4. Edit the “0.10 Chloride Loading (HS)” cycle to utilize the following Cycle Steps and Parameter Values:

Step	Operation	0.2 M Amino Acid Parameter Values	0.5 M Amino Acid Parameter Values
1	Coupling (Chloride Loading)	Microwave Method: Chloride Loading Amino Acid: (from method) Amino Acid Volume: 5 mL Base Bottle Position: Position 21 Base Volume: 2 mL Manifold Wash Volume: 2 mL	Microwave Method: Chloride Loading Amino Acid: (from method) Amino Acid Volume: 2 mL Base Bottle Position: Position 21 Base Volume: 2 mL Manifold Wash Volume: 2 mL
2	Wash	Volume: 4.0 mL Drain Time: 5.0 mL	Volume: 4.0 mL Drain Time: 5.0 mL
3	Wash	Volume: 4.0 mL Drain Time: 5.0 mL	Volume: 4.0 mL Drain Time: 5.0 mL
4	Coupling (Chloride Loading)	Microwave Method: Chloride Loading Amino Acid: (from method) Amino Acid Volume: 5 mL Base Bottle Position: Position 21 Base Volume: 2 mL Manifold Wash Volume: 2 mL	Microwave Method: Chloride Loading Amino Acid: (from method) Amino Acid Volume: 2 mL Base Bottle Position: Position 21 Base Volume: 2 mL Manifold Wash Volume: 2 mL
5	Wash Thru Manifold	Volume: 4.0 mL Drain Time: 5.0 mL	Volume: 4.0 mL Drain Time: 5.0 mL
6	Wash	Volume: 4.0 mL Drain Time: 5.0 mL	Volume: 4.0 mL Drain Time: 5.0 mL
7	Wash	Volume: 4.0 mL Drain Time: 5.0 mL	Volume: 4.0 mL Drain Time: 5.0 mL
8	Wash	Volume: 4.0 mL Drain Time: 5.0 mL	Volume: 4.0 mL Drain Time: 5.0 mL

5. Once all cycle steps and parameters are entered select “Save.”



## Step Two: Apply Cycle to a Liberty Method

1. From the “Edit” tab, select “Liberty methods...” to create a new Liberty Method.
2. Edit the Liberty Method, utilizing the following method options:
  - C-Terminus: Acid
  - Resin-Type: Chloride Loading (or “Trityl DCA Loading” if using software v1.5 or earlier)
  - Resin Cycle: Chloride Loading No Swelling (or “Trityl DCA Loading No Swelling” if using software v1.5 or earlier)
3. In the Amino Acid Cycles grid, double-click the C-terminal amino acid and select the “Chloride Loading” cycle from the drop-down menu, if not already selected.
4. Select “Save” and close the Cycle Editor.
5. Load the Liberty Method into the resin indicator position.
6. To ensure sufficient reagent solutions are prepared, select the “Calculators” tab followed by “Usage Calculator.” See “Step Three: Prepare Reagents” for important reagent preparation tips and considerations.
7. Press play to start the Liberty Method.

## Step Three: Prepare Reagents

1. Weigh Cl-MPA ProTide or Cl-TCP(Cl) ProTide resin and transfer into a clean, dry reaction vessel. Secure the reaction vessel onto the attenuator and place into microwave cavity.

### **WARNING**

The Cl-TCP(Cl) ProTide and Cl-MPA ProTide resins can hydrolyze, reducing yield.

### **NOTE**

If using HT 12/24, load the resin into a clean, dry centrifuge tube and place directly onto the desired HT position. No solvent should be added.

2. Prepare a 1.0 M DIEA + 0.125 M KI solution and load onto Position 21.
  - 2.1. Dissolve anhydrous KI (0.52 g) in DMF (20 mL).
  - 2.2. Add DIEA (4.35 mL) to the KI solution.
  - 2.3. Dilute the solution up to volume of 25 mL.

### **NOTE**

For solution concentration for other synthesis scales, contact: [synthesis.support@cem.com](mailto:synthesis.support@cem.com)

3. If employing Cl-TCP(Cl) ProTide resin, add 0.1 equivalents of DIEA to the HOBt or Oxyma solution.

### **WARNING**

Cl-TCP(Cl) ProTide resin is hyper-acid sensitive. Addition of DIEA to the HOBt or Oxyma solution is necessary for prevention of premature cleavage, increasing synthesis yield.

4. Prepare any additional necessary reagents, load onto instrument, and begin the Liberty Method.

## Optional: Protected Cleavage with CI-TCP(CI) ProTide

Resins with hyper-acid sensitive linkers like CI-TCP(CI) ProTide can be treated with 1% TFA to cleave the peptide from resin while the protecting groups stay intact.

### **i** NOTE

The Liberty Blue is not equipped to perform peptide cleavage; all cleavage procedures must be performed off-system.

The procedure below is for peptide cleavage at the 0.10 mmol scale. Adjust as necessary.

1. Transfer the peptidyl resin from the Liberty Blue reaction vessel to a filtered syringe vessel.
2. Rinse the resin thoroughly with DCM.
3. Add 1% TFA in DCM (double the resin bed volume) to the syringe vessel.
4. After 2 min at room temperature, collect the filtrate in a centrifuge tube containing a solution of 10% pyridine in MeOH (3 mL).
5. Repeat steps 3 and 4 four additional times.
6. Combine all filtrates, transferring quantitatively with MeOH and concentrate the solution (< 3 mL).
7. Transfer the solution to a clean 50-mL centrifuge tube with a small amount of MeOH (< 2 mL).
8. Dilute the solution to 35 mL with ice-cold deionized water to precipitate the protected peptide.
9. Centrifuge the peptide solution for 5 min at 3500 rpm, or until completion of peptidyl pellet formation. Upon decantation of the supernatant the peptide is ready for lyophilization and/or analysis.