



Procedure

Weigh 350 mg of the drillings into the digestion vessel. Add 6 mL of HCl, 2 mL of HNO₃ and 2 mL of HF. Gently swirl the mixture and wait approximately 15 minutes before closing the vessel. Allow any initial reaction to subside before sealing vessel.

Notes

This application can only be run in the iPrep vessel.

HF should be added slowly and carefully to the sample.

This procedure uses hydrofluoric acid. If it is necessary to complex the residual hydrofluoric acid or redissolve insoluble fluorides formed, an additional complexation step with boric acid or other neutralization step should be used.

Recommended Equipment

MARS 6 iWave

Recommended Vessels

iPrep

Reagents

HCl
HNO₃
HF

Max Sample Weight

350 mg

Sample Type

Organic

Control Type

Ramp to Temperature

Method Type

One Touch

Heating Program

Stage	Temp (°C)	*Ramp (mm:ss)	Hold (mm:ss)	Pressure (psi)	* Power (W)	Stirring
1	200	30:00	15:00	N/A	700-1800	Off

* Ramp times and power may vary depending on the type and number of vessels.

Results

Sample was clear, emerald green in color with no visible particles upon dilution to 50 mL.

General Precaution

- This procedure is a reference point for sample digestion using a CEM system and may need to be modified or changed to obtain the required results on your sample.
- If using a vessel other than the recommended choice, adjust sample size and pressure limit to values appropriate for the vessel chosen.
- Manual venting of CEM vessels should be performed when wearing hand/eye/body protection and when the vessel contents are at or below room temperature to avoid the potential for chemical burns. Always point the vent hole away from the operator.
- If programming as One Touch, the ramp time and power will be automatically determined based on the number and type of vessels detected.
- If using HF, follow restrictions listed in HF Addendum.