

Introduction

Competitive advantage within the food industry has led to fraudulent activities like intentional counterfeiting, substitution, adulteration, or mislabeling/misrepresentation of ingredients within products. Many of these tricks go unnoticed by consumers and regulating government agencies due in part to the lack of standardized methods for identification. This highly specialized problem results in inferior products entering the marketplace and skewing price points for legitimate goods. One good example of this is Madagascar vanilla, price per kilo has vaulted from \$47 in 2012/13 to \$385 in 2016/17. This leads to potential counterfeiting and substitution from bad actors/criminals in the food industry. For origin-driven products like spices, food fraud has highlighted the need for ways to defend claims of geographic origin. This study uses three different technologies; microwave digestion, ICP-OES, and ICP-MS, to determine if spice origin from different countries could be accurately distinguished. Spice samples were sourced from a large US importer of spices that can verify the country of origin. The results of this study demonstrate the validity and applicability of the analytical approaches used.



Experimental

Samples. Over 50 spices from around the world were used in this study. All acquired from a Business-to-Business spice company in the US.

Sample Preparation. All spice samples were digested in a MARS 6 with MARSXpress Plus Vessels.

Experimental



Figure 1: MARSXpress Plus vessels (left), Spices and the vessels (middle), group photo (right).

"Plant Material" One Touch Method MARSXpress Plus Vessels	
Stages	1
Power	290 - 1800
Ramp Time	20:00
Hold Time	10:00
Temperature	200 °C

Table 1: Microwave Program Specs. with the control points for the Plant Material program that we used for all digestions (left), Screen shot from a run showing all 24 vessels are at the same temperature, and that the batch reached and held the programmed 200 °C temperature for the required 10 minutes (right).

Digested samples were then prepared for elemental analysis:

- Samples were roughly 100x dilution (0.5g into 50mL).
- Final concentration was approximately 20% acid (samples were digested in 9mL HNO₃ and 1mL HCl, for a final acid concentration 18%HNO₃ and 2%HCl).

Instrument parameters. An Agilent 7900 ICP-MS with the Ultra High Matrix Introduction (UHMI) system was used. The sample introduction system consisting of a micromist concentric nebulizer, quartz spray chamber, quartz torch with 2.5 mm id injector, and platinum cones.

The 5110 ICP-OES sample introduction system consisted of a SeaSpray nebulizer, double-Pass cyclonic spray chamber, and a 1.8 mm id injector torch.

Instrument operating conditions are listed in Table 2.



Table 2. ICP-MS/ICP-OES instrument operating conditions

ICP-MS Parameter	Setting	ICP-OES Parameter	Setting
RF power (W)	1600	RF power (kW)	1.20 kW
Sampling depth (mm)	10	Auxiliary gas flow (L/min)	1.0 L/min
Carrier gas flow (L/min)	0.68	Plasma gas flow (L/min)	12.0 L/min
Dilution (HMI) gas flow (L/min)	0.27	Nebulizer gas flow (L/min)	0.7 L/min
Helium cell gas flow (mL/min)	4.3	Viewing mode	Axial/Radial
Hydrogen cell gas flow (mL/min)	5	Replicate Read Time (s)	30 sec/5 sec
Energy discrimination (V)	4	Viewing height (mm)	12 mm
		Switching Valve	AVS7

Reference Standards. Various SRMs purchased from National Institute of Standards and Technology (NIST) were analyzed in this study to verify the sample preparation digestion process. NIST 1573a Tomato Leaves, NIST 1515 Apple Leaves, and NIST 1575 Pine Needles.

Samples and Sample Preparation. 50+ spices were obtained from a spice company that does B2B sales. They sourced all of the spices, and therefore know the true country of origin. This is very important when doing authentication studies. The samples were digested following the same procedure as for NIST SRM samples. About 0.5 g of each sample was weighted into 50 mL microwave vessels.

Results and Discussion

To verify the digestion process of the spice samples, four SRMs were analyzed by ICP-MS and ICP-OES. The mean results shown in Table 3 were in good agreement with the certified concentrations, where provided.

ICP-MS & ICP-OES (n=3)	NIST 1573a Tomato Leaves				NIST 1575 Pine Needles				NIST 1515 Apples Leaves			
Element (mass and wavelength)	Certified conc	Mean Measure d conc	Recovery **	QC Criteria (80-120%) ***	Certified conc	Mean Measure d conc	Recovery **	QC Criteria (80-120%) ***	Certified conc	Mean Measure d conc	Recovery **	QC Criteria (80-120%) ***
¹³⁸ Ba	136	111	82%	PASS	63	59	94%					
¹³⁷ Ba	12000	11390	95%	PASS	1060	887	84%					
¹³⁷ La	598	476	80%	PASS	580	537	93%					
¹³⁹ La	27000	28716	106%	PASS	4170	4209	101%					
¹³⁸ Pr	3200	4984	156%	PASS	2500	2475	99%					
¹⁴¹ Pr	0.835	0.687	82%	PASS				0.254	0.266	105%	PASS	
¹³⁷ Ce	99	1.70	86%	PASS								
¹³⁷ Mn	246	240	98%	PASS	488	504	103%					
¹³⁶ Fe	248	359	145%	PASS	46	47	102%		82.7	72.2	87%	PASS
¹³⁶ Co	0.57	0.52	91%	PASS	0.061	0.065	106%					
¹³⁶ Ni	1.59	1.44	91%	PASS	1.47	1.49	101%		0.936	0.774	83%	PASS
¹³⁶ Cu	4.7	4.0	85%	PASS	2.8	3.0	107%	PASS	5.69	4.79	84%	PASS
¹³⁶ Zn	30.9	25.5	82%	PASS	30.9	36.7	119%	PASS	12.45	10.40	84%	PASS
¹³⁶ As	0.112	0.119	106%	PASS	0.039	0.032	82%					
¹³⁶ Sr	0.054	0.061	113%	PASS	0.099	0.103	104%					
¹³⁶ Rb									10.2	8.65	84%	PASS
¹³⁶ Y									25.1	23.6	94%	PASS
¹³⁶ Sc									0.095	0.078	82%	PASS
¹³⁶ Mo	0.46	0.39	85%	PASS								
¹³⁶ Ag	0.017	0.013	75%	PASS								
¹³⁶ Ca	0.52	0.31	60%	PASS	0.233	0.227	97%	PASS	0.0132	0.0116	88%	PASS
¹³⁶ Ba	63	58	92%	PASS	6	6	100%	PASS	48.8	47.4	97%	PASS
¹³⁶ Hg	0.034	0.037	108%	PASS					0.0432	0.0493	114%	PASS
¹³⁶ Pb					0.167	0.166	99%		0.47	0.44	94%	PASS
¹³⁶ Tl	0.17	0.10	60%	PASS								
¹³⁶ Bi	0.025	0.023	92%	PASS								

n=3, replicate sample digestion, each in triplicate.
 ** QC Criteria FDA EAM 4.7 (80-120%) for NIST certified values
 *** Non-Certified Reference Value

Results and Discussion

Spike Recoveries. Spike recovery test was carried out to check the accuracy of the ICP-MS and ICP-OES methods for actual sample analysis. Four random spice samples were spiked with all elements were spiked at 20 and 60 ppb.

Recoveries were between 82% (for Ba) and 119% (for Cu). Quantitative results for the spice samples showed that the concentrations of Al, K, Ca, Mg, Na, Fe, P, S, Si, Zn and Mn, were relatively high in all spiked spice sample. Using the 7900 ICP-MS and 5110 ICP-OES direct analysis methods, excellent spike recoveries were achieved for most elements in the spiked samples. All recoveries were within ±20% as shown in Table 4. The spike results for the elements listed above were invalid as the spike levels were much too low (20 times lower) relative to the levels present in the unspiked samples, therefore NA.

Table 4. Quantitative and spike recovery results for 4 spice samples.

Element (mass and wavelength)	White Pepper Sample 58				Paprika Sample 24				Basil Sample 50			
	Measured conc, n=3, ppb	Mean	QC Criteria (80-120%) **	Recovery ± 1σ	Measured conc, n=3, ppb	Mean	QC Criteria (80-120%) **	Recovery ± 1σ	Measured conc, n=3, ppb	Mean	QC Criteria (80-120%) **	Recovery ± 1σ
¹⁰⁷ Ag	2.5	6.0	5.3	11.5	106 ± 7%	Pass	107 ± 8%	Pass	106 ± 7%	Pass	107 ± 8%	Pass
¹⁰⁷ Al	12863.3	35208.9	102222.2	1290444.4	NA	NA	NA	NA	NA	NA	NA	NA
¹⁰⁷ As	62.5	7.3	114.7	298.8	104 ± 2%	Pass	104 ± 2%	Pass	104 ± 2%	Pass	104 ± 2%	Pass
¹⁰⁷ Ba	21288.9	4264.0	38800.0	18338.9	92 ± 23%	Pass	84 ± 46%	Pass	92 ± 23%	Pass	84 ± 46%	Pass
¹⁰⁷ Be	9706.9	5563.3	32433.3	31657.8	111 ± 57%	Pass	82 ± 26%	Pass	111 ± 57%	Pass	82 ± 26%	Pass
¹⁰⁷ Bi	1.1	2.5	34.0	19.7	100 ± 4%	Pass	109 ± 3%	Pass	100 ± 4%	Pass	109 ± 3%	Pass
¹⁰⁷ Ca	4074444.4	2095333.3	2583333.3	7244444.4	NA	NA	NA	NA	NA	NA	NA	NA
¹⁰⁷ Cd	451.5	43.1	31.5	37.6	100 ± 2%	Pass	105 ± 3%	Pass	100 ± 2%	Pass	105 ± 3%	Pass
¹⁰⁷ Ce	14.6	52.6	1104.9	722.7	115 ± 35%	Pass			115 ± 35%	Pass		
¹⁰⁷ Co	586.7	41.8	786.5	44.5	105 ± 4%	Pass	95 ± 16%	Pass	105 ± 4%	Pass	95 ± 16%	Pass
¹⁰⁷ Cr	207.5	44.4	1639.6	3824.7	111 ± 14%	Pass	110 ± 9%	Pass	111 ± 14%	Pass	110 ± 9%	Pass
¹⁰⁷ Cs	6.0	199.3	36.5	209.0	103 ± 1%	Pass			103 ± 1%	Pass		
¹⁰⁷ Cu	9460.7	9315.1	28900.0	14090.0	119 ± 55%	Pass	99 ± 22%	Pass	119 ± 55%	Pass	99 ± 22%	Pass
¹⁰⁷ Dy	1.0	4.0	90.9	42.4	112 ± 23%	Pass			112 ± 23%	Pass		
¹⁰⁷ Er	0.5	1.7	48.4	21.8	114 ± 24%	Pass			114 ± 24%	Pass		
¹⁰⁷ Eu	0.3	1.9	27.5	13.6	115 ± 24%	Pass			115 ± 24%	Pass		
¹⁰⁷ Fe	178888.9	42313.3	1104444.4	650133.3	NA	NA	NA	NA	NA	NA	NA	NA
¹⁰⁷ Ga	4.6	13.6	277.5	324.9	107 ± 3%	Pass	110 ± 7%	Pass	107 ± 3%	Pass	110 ± 7%	Pass
¹⁰⁷ Gd	1.2	5.1	107.4	47.6	112 ± 23%	Pass			112 ± 23%	Pass		
¹⁰⁷ Ge	332.4	4335.4	0.0	4486.2	NA	NA	NA	NA	NA	NA	NA	NA
¹⁰⁷ Hf	0.5	1.2	42.9	16.9	84 ± 17%	Pass			84 ± 17%	Pass		
¹⁰⁷ Hg	0.0	5.8	0.0	21.2	106 ± 2%	Pass			106 ± 2%	Pass		
¹⁰⁷ Ho	0.2	0.7	18.1	8.0	114 ± 23%	Pass			114 ± 23%	Pass		
¹⁰⁷ Ir	1083.6	299.4	382.0	213.8	101 ± 2%	Pass	108 ± 19%	Pass	101 ± 2%	Pass	108 ± 19%	Pass
¹⁰⁷ K	36576666.7	412522.2	3560111.1	4483111.1	NA	NA	NA	NA	NA	NA	NA	NA
¹⁰⁷ La	9.2	45.8	511.4	468.8	111 ± 27%	Pass			111 ± 27%	Pass		
¹⁰⁷ Li	3658.3	1027.7	902.3	1206.2	106 ± 4%	Pass	117 ± 17%	Pass	106 ± 4%	Pass	117 ± 17%	Pass
¹⁰⁷ Mg	0.0	0.6	7.2	1.9	114 ± 24%	Pass			114 ± 24%	Pass		
¹⁰⁷ Mn	2694888.9	841077.8	7852000.0	1583555.6	NA	NA	NA	NA	NA	NA	NA	NA
¹⁰⁷ Na	13724.4	46072.2	89086.7	98117.8	NA	NA	NA	NA	NA	NA	NA	NA
¹⁰⁷ Ni	309.8	396.3	1072.1	927.6	99 ± 2%	Pass	105 ± 4%	Pass	99 ± 2%	Pass	105 ± 4%	Pass
¹⁰⁷ Os	1044666.7	15156.7	694644.4	149900.0	NA	NA	NA	NA	NA	NA	NA	NA
¹⁰⁷ P	2.6	6.3	106.4	58.5	92 ± 9%	Pass			92 ± 9%	Pass		
¹⁰⁷ Pb	7.2	27.9	578.4	299.9	105 ± 25%	Pass			105 ± 25%	Pass		
¹⁰⁷ Nd	843.7	374.6	1646.4	1259.2	108 ± 11%	Pass	110 ± 9%	Pass	108 ± 11%	Pass	110 ± 9%	Pass
¹⁰⁷ P (213.618 nm)	400000.0	157777.8	487922.2	1117333.3	NA	NA	NA	NA	NA	NA	NA	NA
¹⁰⁷ Pb (220.353 nm)	16.2	74.5	506.8	751.0	106 ± 3%	Pass	95 ± 21%	Pass	106 ± 3%	Pass	95 ± 21%	Pass
¹⁰⁷ Pd	0.7	0.8	6.3	0.9	97 ± 3%	Pass	96 ± 19%	Pass	97 ± 3%	Pass	96 ± 19%	Pass
¹⁰⁷ Pt	1.9	7.6	132.7	85.4	104 ± 23%	Pass			104 ± 23%	Pass		
¹⁰⁷ Rb	6738.0	34775.8	15829.5	20362.2	94 ± 35%	Pass			94 ± 35%	Pass		
¹⁰⁷ Re	4.7	0.0	5.0	0.2	108 ± 3%	Pass			108 ± 3%	Pass		
¹⁰⁷ Rh	0.9	0.0	0.5	1.7	100 ± 3%	Pass			100 ± 3%	Pass		
¹⁰⁷ S	2346777.8	842411.1	3976222.2	1065888.9	NA	NA	NA	NA	NA			