

Microwave Assisted SPPS of Hindered, Non-Standard Amino Acids




Summary


- Microwave-enhanced SPPS enables conventionally-difficult couplings of bulky amino acids, like Aib and *N*-Me-A, to occur quickly and efficiently.
- Syntheses of acyl carrier protein derivatives VQ(Aib)₂IDYING-OH and VQ(*N*-Me-A)₂(*N*-Me-A)IDYING-OH are completed in under 2 h and in 95% and 86% purities respectively
- Synthesis of GEQKLG(Aib)₂AibASEESLG-NH₂ is completed in under 3 hr with an 89% purity.

Introduction

Hindered, non-standard amino acids such as α -aminoisobutyric acid (Aib) and *N*-methyl alanine ((*N*-Me)-A) (**Figure 1**) can be found in many biologically relevant compounds.¹⁻³ The synthesis of peptides including Aib or *N*-methylated amino acids has proved challenging, however; the steric hindrance introduced by the second methyl group, whether on the α -carbon or the amide nitrogen, makes coupling these amino acid derivatives difficult in conventional SPPS.



α -Aminoisobutyric Acid



N-Methyl Alanine

Figure 1 Sterically-Hindered, Non-Standard Amino Acids

Through the use of microwave-enhanced SPPS, though, difficulties associated with hindered, non-standard amino acids have been minimized. The employment of microwave energy in SPPS drives conventionally-difficult couplings of bulky amino acids, like Aib and *N*-methyl alanine, quickly and efficiently to completion.^{4,5}

Materials and Methods

Reagents

N- α -Fmoc- α -aminoisobutyric acid was obtained from AnaSpec (Freemont, CA). Fmoc-*N*-Me-Ala-OH was obtained from Peptides International (Louisville, KY). All other amino acids were obtained from CEM Corporation (Matthews, NC) and contained the following side chain protecting groups: Asn(Trt), Asp(OMpe), Gln(Trt), Glu(OtBu), Lys(Boc), Ser(OtBu), and Tyr(tBu). Oxyma Pure and Rink Amide ProTide™ LL resin were obtained from CEM Corporation (Matthews, NC). *N,N*-Diisopropylcarbodiimide (DIC) was obtained from CreoSalus (Louisville, KY). Fmoc-Gly-Wang Resin LL was obtained from NovaBiochem (St. Louis, MO). Piperidine was obtained from Alfa Aesar (Ward Hill, MA). Trifluoroacetic acid (TFA), 3,6-dioxa-1,8-octanedithiol (DODT), triisopropylsilane (TIS), and acetic acid were obtained from Sigma-Aldrich (St. Louis, MO). Dichloromethane (DCM), *N,N*-dimethylformamide (DMF), and anhydrous diethyl ether (Et₂O) were obtained from VWR (West Chester, PA). HPLC-grade water (H₂O), and HPLC-grade acetonitrile (MeCN) were obtained from Fisher Scientific (Waltham, MA).

Peptide Synthesis: GEQKLG(Aib)₂AibASEESLG-NH₂

The peptide was prepared at 0.1 mmol scale using the CEM Liberty Blue™ automated microwave peptide synthesizer

on Rink Amide ProTide LL resin (0.18 meq/g substitution). Deprotection was performed with piperidine and Oxyma Pure in DMF. Coupling reactions were performed with DIC in DMF, Oxyma Pure in DMF, and a 5-fold excess of Fmoc-AA-OH. Cleavage was performed using the CEM Razor high-throughput peptide cleavage system with TFA/H₂O/TIS/ DODT. Following cleavage, the peptide was precipitated in Et₂O and lyophilized overnight.

Peptide Synthesis: VQAibAibIDYING-OH

The peptide was prepared at 0.1 mmol scale using the CEM Liberty Blue automated microwave peptide synthesizer on Fmoc-Gly-Wang LL Resin (0.33 meq/g substitution). Deprotection was performed with piperidine and Oxyma Pure in DMF. Coupling reactions were performed with DIC in DMF, Oxyma Pure in DMF, and a 5-fold excess of Fmoc-AA-OH. Cleavage was performed using the CEM Razor high-throughput peptide cleavage system with TFA/H₂O/TIS/ DODT. Following cleavage, the peptide was precipitated in Et₂O and lyophilized overnight.

Peptide Synthesis: VQ(N-Me-A)(N-Me-A)IDYING-OH

The peptide was prepared at 0.1 mmol scale using the CEM Liberty Blue automated microwave peptide synthesizer on Fmoc-Gly-Wang LL Resin (0.19 meq/g substitution). Deprotection was performed with piperidine and Oxyma Pure in DMF. Coupling reactions were performed with DIC in DMF, Oxyma Pure in DMF, and a 5-fold excess of Fmoc-AA-OH. Cleavage was performed using the CEM Razor high-throughput peptide cleavage system with TFA/H₂O/TIS/ DODT. Following cleavage, the peptide was precipitated in Et₂O and lyophilized overnight.

Peptide Analysis

The peptides were analyzed on a Waters Acquity UPLC system with PDA detector equipped with an Acquity UPLC BEH C8 column (1.7 mm and 2.1 x 100 mm). The UPLC system was connected to a Waters 3100 Single Quad MS for structural determination. Peak analysis was achieved on Waters MassLynx software. Separations were performed with a gradient elution of 0.1% TFA in (i) H₂O and (ii) MeCN.

Results

Microwave-enhanced SPPS of GEQKLGAIbAibAibASEEDLG-NH₂ on the Liberty Blue automated microwave peptide synthesizer produced the target peptide in 89% purity (**Figure 2**).




Figure 2 UPLC Chromatogram of GEQKLGAIbAibAibASEEDLG-NH₂

Microwave-enhanced SPPS of VQAibAibIDYING-OH on the Liberty Blue automated microwave peptide synthesizer produced the target peptide in 95% purity (**Figure 3**).




Figure 3 UPLC Chromatogram of VQAibAibIDYING-OH

Microwave-enhanced SPPS of VQ(N-Me-A)(N-Me-A)IDYING-OH on the Liberty Blue automated microwave peptide synthesizer produced the target peptide in 86% purity (**Figure 4**).




Figure 4 UPLC Chromatogram of VQ(N-Me-A)(N-Me-A)IDYING-OH

Conclusion

Microwave-enhanced SPPS enables conventionally-difficult couplings of bulky amino acids, like Aib and *N*-Me-A, to occur quickly and efficiently. Though conventional synthesis produces GEQKLG Aib Aib Aib ASEEDLG-NH₂ in 40 h and < 10% purity, microwave-enhanced SPPS produces the target peptide in under 3 h and in 89% purity. Additionally, syntheses of acyl carrier protein derivatives VQAib Aib IDYING-OH and VQ(*N*-Me-A)(*N*-Me-A)IDYING-OH are completed in under 2 h and in 95% and 86% purities respectively. Microwave-enhanced SPPS has proven an effective tool in minimizing the difficulties associated with hindered, non-standard amino acids in SPPS.

References

- (1) Mueller, P.; Rudin, D. O. *Nature*. **1968**, *217*, 713–719.
- (2) Rebuffat, S.; Goulard, C.; Hlimi, S.; Bodo, B. *J. Pept. Sci.* **2000**, *6*, 519–533.
- (3) Ahmed, G.; Elger, W.; Meece, F.; Nair, H. B.; Schneider, B.; Wyrwa, R.; Nickisch, K. *Bioorg. Med. Chem.* **2017**, *25*, 5569–5575.
- (4) Collins, J. M. Microwave-Enhanced Synthesis of Peptides, Proteins, and Peptidomimetics. In *Microwaves in Organic Synthesis*, Third Edition; de la Hoz, A.; Loupy, A.; Wiley-VCH: Weinheim, 2012; 897–959.
- (5) Ben Haj Salah, K.; Inguibert, N. *Org. Lett.* **2014**, *16* (6), 1783–1785.

United States (Headquarters)

800-726-3331
704-821-7015
Fax: 704-821-7894
info@cem.com

France

33 (01) 69 35 57 80
Fax: 33 (01) 60 19 64 91
info.fr@cem.com

Germany, Austria, Switzerland

(49) 2842-9644-0
Fax: (49) 2842-9644-11
info@cem.de

Ireland

+353 (0) 1 885 1752
Fax: +353 (0) 1 885 1601
info.ireland@cem.com

Italy

(39) 35-896224
Fax: (39) 35-891661
info.srl@cem.com

Japan

+81-3-5793-8542
Fax: +81-3-5793-8543
info@cemjapan.co.jp

United Kingdom

(44) 1280-822873
Fax: (44) 1280-822873
info.uk@cem.com

www.cem.com

© 2018 CEM Corporation
All rights reserved. This may not be reproduced or published without written permission from CEM.