


Esterification Experiment

Introduction


Esters are organic compounds derived from carboxylic acids in which the acidic hydrogen has been replaced with an organic group. They are known for their pleasant odors and comprise the major flavor components of a number of fruits. In addition, they play a significant role in primary metabolism and other biochemical processes.

Esters can be synthesized in a number of different ways. One common method is the acid catalyzed condensation of a carboxylic acid and an alcohol. The mechanism involves an initial protonation of the carboxylic acid to form a resonance-stabilized cation. The addition of the alcohol to the cationic carbon and loss of one of the protons on the oxygen of the alcohol form a species called a geminal diol. Protonation of one of the hydroxyl groups in this diol, followed by loss of water, yields the ester product. To obtain significant yields of the product, the equilibrium must be shifted by either added excess of one of the starting materials or by removing water and/or the ester product.



In this experiment, benzyl alcohol is used in the esterification of ethanoic (acetic) acid to give a peach-smelling ester. Using conventional heating, the esterification reaction reaches equilibrium after a few hours of reflux. In this experiment, microwave heating in the CEM Discover SP™ microwave system is used to significantly accelerate the reaction. Both excess acetic acid and the removal of the water are utilized to drive the reaction to completion.

Materials and Methods



Reagent	MW (g/mol)	mmol	Mass (g)	Density (g/mL)	Volume (mL)
Glacial acetic acid	60	35	2.1	1.049	2.0
Benzyl Alcohol	108	12	1.3	1.045	1.2
Sulfuric Acid					10 drops
Silica beads			0.20		

Procedure for benzyl ethanoate (benzylacetate)

Glacial acetic acid (2.0 mL, 35 mmol), benzyl alcohol (1.2 mL, 12 mmol, 1 eq), concentrated sulfuric acid (10 drops), and silica beads (0.20 g) are placed in a 10-mL Discover SP™ reaction vessel containing a stir bar. The reaction vessel is

sealed with a cap and then placed in the Discover SP™ microwave cavity. **Program a Dynamic Method to heat the reaction mixture to 130 °C and hold for 5 minutes.** After the reaction had completed and the vessel has cooled to below 50 °C, remove the vessel from the microwave.

Sodium bicarbonate (NaHCO_3) (10 mL of a 10% solution) is placed in a separatory funnel. The reaction mixture is transferred from the microwave vessel into the funnel with a pipette. This will lead to a very effervescent reaction as the excess acetic acid is neutralized. Water is added at this point until any visible solid is dissolved. Diethyl ether (5 mL) is added to the organic layer and the solution swirled. The aqueous layer is then removed. An additional 5 mL of sodium bicarbonate is added and the solution swirled again. After the layers have separated, the organic layer is removed to a small Erlenmeyer flask. The aqueous layer is extracted two more times with 5 mL of diethyl ether, each time. Remove the organic layer after each extraction and combine it with the first ether layer. The combined organic layers are returned to the empty separatory funnel and are washed with saturated NaCl solution (5 mL). The organic layer is placed in a clean Erlenmeyer flask and then dried over sodium sulfate (Na_2SO_4). The dried organics are transferred into a pre-weighed flask. The diethyl ether is removed on a rotary evaporator or under a stream of nitrogen. The flask plus product is weighed and the yield calculated for the reaction. The refractive index is determined and compared to that in the literature. The chemical structure index is determined and compared to that in the literature. The chemical structure of the product is confirmed by IR, ^1H -NMR and/or ^{13}C -NMR spectroscopy.

"Clean, Fast Organic Chemistry: Microwave-assisted Laboratory Experiments." Nicholas Leadbeater, PhD. and Cynthia McGowan, Ph.D.



United States (Headquarters):

CEM Corporation
PO Box 200
Matthews, NC 28106
800-726-3331
704-821-7015
Fax: 704-821-7894
info@cem.com
cem.com

France:

CEM μ Wave S.A.S.
Immeuble Ariane
Domaine Technologique de Saclay
4, rue Rene' Razel, 91892 ORSAY Cedex
33 (01) 69 35 57 80
Fax: 33 (01) 60 19 64 91
info.fr@cem.com

Germany, Austria, Switzerland:

CEM GmbH
Carl-Friedrich-Gauss-Str.9,
47475 Kamp-Lintfort
(49) 2842-9644-0
Fax: (49) 2842-9644-11
info@cem.de

Ireland:

CEM Technology (Ireland) Ltd.
Sky Business Centre,
9a Plato Business Park,
Damastown, Dublin 15
+353 (0) 1 885 1752
Fax: +353 (0) 1 885 1601
info.ireland@cem.com

Italy:

CEM S.R.L.
Via Dell' Artigianato, 6/8
24055 Cologno al Serio (Bg)
(39) 35-896224
Fax: (39) 35-891661
info.srl@cem.com

Japan:

CEM Japan K.K.
2-18-10 Takanawa, Minato-ku, Tokyo
108-0074
+81-3-5793-8542
Fax: +81-3-5793-8543
info@cemjapan.co.jp

United Kingdom:

CEM Microwave Technology Ltd.
2 Middle Slade, Buckingham Industrial
Estate, Buckingham MK181WA
(44) 1280-822873
Fax: (44) 1280-822873
info.uk@cem.com