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**Application Note**  
The Extraction of Antioxidants from Food Packaging Film

# The Extraction of Antioxidants from Food Packaging Film

## Abstract

Polymer materials are an integral part of every day life, especially food packaging. Often, antioxidants are added to food packaging materials to help aid in their stability. However, these antioxidants can leach out of the materials into food, demonstrating the importance of being aware of the levels of antioxidants in the polymer materials in food packaging. In this work, Irganox 1076 and Irganox 1010 were extracted from food film provided by a local manufacturer, using three CEM extraction systems, the MARS 6™ microwave system, the Discover Prep™ microwave system, and the EDGE® automated extraction system. The MARS 6 method has been long proven in the industry and is known to deliver accurate results. The results obtained for the extractions completed by Discover Prep and the EDGE were compared to the MARS 6 results, and it was found that the Discover Prep and the EDGE extracted the additives from food

## Materials and Methods

### Reagents

Food film was provided by a local manufacturer. ASTM D6042-96 Calibration Mix was purchased from Restek. ACS-grade 2-propanol and HPLC-grade 2-propanol were purchased from Sigma. 5-mL syringes were purchased from Fisher, and 0.2-micron PTFE Luer lock syringe filters were purchased from Phenomenex. Other materials were provided by CEM.

### MARS Extraction

A portion of 1 g of food film, cut into 0.5-cm wide squares, was weighed directly into a 110 mL MARSXpress™ Plus vessel, and 20 mL of ACS-grade 2-propanol were added. Samples were prepared in quadruplicate. Each vessel was loaded into

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[cem.com](http://cem.com)

United States Headquarters: +1 (704) 821-7015 | [info@cem.com](mailto:info@cem.com)  
For distributors and subsidiaries in other regions, visit [cem.com/contact](http://cem.com/contact)

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