



Discover[®] SPPS

Microwave Synthesizer





WASTE

Discover SPPS Solid Phase Peptide Synthesis System

Advantages of Microwave Heating

Benchtop solid phase peptide synthesis (SPPS) is a time intensive process that often leads to challenges in synthesis and purification. Microwave heating can reduce reaction times and improve purity in peptide synthesis.



Accelerated Reaction Rates

Improved Synthetic Purity



Small Benchtop Footprint



Economically Priced

Pioneers in Microwave SPPS

CEM developed and continues to innovate in microwave assisted SPPS. Microwave heating provides rapid, responsive, and higher efficiency heating that leads to optimized methodology for high speed and purity. The Discover SPPS leverages this legacy for high level performance in a manual synthesizer.



SPPS Method Development

The Discover SPPS empowers peptide scientists to investigate each step of a peptide synthesis through easy access to the reaction mixture after every reaction. The result of each chemical step can be aliquoted or paused for analysis. The Discover SPPS can be used to test reaction conditions before committing to a full peptide synthesis run.

Meeting All Your Needs



Fast, Flexible, High-Purity Synthesis

	Liberty Lite 2.0	Liberty Blue 2.0
Cycle Time (at 0.1mmol)	12 minutes	4 minutes
Waste/Cycle (at 0.1mmol)	32 mL	16 mL
Scale Range	0.005 - 5 mmol	0.005 - 5 mmol
Ultra Purity with HeadSpace Flushing	No	Yes
Amino Acid Positions	27	27
Multiple Peptide Synthesis	Up to 4 (HT4)	Up to 4 (HT4) Up to 12 (HT12)

Upgrade to Automation

Expand to our Liberty 2.0 Series upgrades for unprecedented speed and efficiency. Access full cycle automation and high-throughput with 4, 12, and 24 peptides.

.....

Hands on Teaching

Engage students with SPPS by taking them through the process step by step. Utilize rapid microwave synthesis times to streamline time spent in lab. Accelerate coupling and deprotections to illustrate important chemistry and refinement in the peptide synthesis process.



Specifications

Vessel volume: Chemistry: Depro + Coupling Time: Temperature Sensor: 25 mL, optimized for 0.1 mmol Fmoc, t-Boc 8 minutes Fiber Optic Maximum Temperature: Reagent Transfer: Cleavage Time: Upgrades: 75 °C Manual or Vacuum 30 minutes MAOS, Liberty Blue 2.0, Liberty Prime 2.0





We Simplify Science

cem.com



United States (Headquarters): 800-726-3331 | info@cem.com For distributors and subsidiaries in other regions, visit cem.com/contact