

# The Elimination of PFAS Background on the EDGE through Material Substitutions



## Abstract

Per- and polyfluoroalkyl substances (PFAS) are a group of compounds that are used in many different industries because of their favorable properties, such as heat and grease resistance. Their resistance and prevalence in manufacturing has presented hurdles for their detection, as they are used in the manufacture of materials typically found in LC-MS systems. This has caused LC-MS companies to suggest the substitution of these materials. Similarly, the EDGE<sup>®</sup>, an automated extraction system, was modified to prevent the contamination of the resultant extracts by PFAS. The cleanliness of the EDGE itself after a series of washes and the background of the S1 Q-Disc<sup>®</sup> stack (C9+G1+C9 sandwich), Q-Matrix Hydra<sup>™</sup>, and the Q-Cup<sup>®</sup> used by the EDGE were examined. It was found that the substitutions provided a PFAS-free background on the EDGE, and the Q-Discs, Q-Matrix Hydra, and Q-Cups were PFAS-free. Thus, the EDGE is an excellent choice for laboratories that would like to automate their PFAS extractions without subjecting their samples to unwanted contamination.

## Introduction

Per- and polyfluoroalkyl substances (PFAS) are a class of man-made chemicals that are used in a wide variety of industries. They are resistant to high temperatures, stains, and grease, allowing for the performance of products such as nonstick cookware, firefighting foam, and stain-resistant fabrics and furnishings. PFAS molecules are made of linked carbon atoms with fluorine atoms branching off of the main chain. The strong bond of carbon and fluorine makes these compounds extremely stable. Their stability, and widespread use, presents difficulties

for their analytical measurement, as they often leach out of the common materials used in LC-MS systems. Because of this, LC-MS companies recommend substituting these components with other materials, such as replacing Teflon-coated solvent lines with PEEK tubing, to remove the background contamination.

The EDGE is an automated extraction system that has been shown to rapidly and efficiently extract PFAS molecules from samples. To ensure that the EDGE was free of any PFAS contamination, substitutions for potentially contaminated materials were made, extra precautions were taken to prevent contamination, and consumables and Q-Cups were assessed for their contamination of PFAS content. The system was found to be free of any PFAS contamination, and the consumables and Q-Cups were also found to have no contamination. Thus, the modifications done to the EDGE ensure a clean background for the extraction of PFAS from a variety of matrices.

## Materials and Methods

### EDGE Modifications

Typically, the EDGE has solvent tubing lines, a dispense line, and other lines along its fluidic pathway made of polytetrafluoroethylene (PTFE). However, PFAS compounds, such as perfluorooctanoic acid, are involved in the manufacture of PTFE. This warrants the replacement of these lines. PEEK tubing was substituted for the PTFE solvent tubing, as the PEEK tubing is flexible enough to be used for this purpose. The PTFE tubing used in the dispense line and other lines along the

pathway were substituted with polypropylene tubing because these are moving components that require more flexibility than PEEK tubing allows. Solvent caps containing PTFE seals were also substituted with polypropylene caps.

As a further safeguard against contamination, extra care was taken in a variety of ways. A side enclosure was utilized on the EDGE to prevent air contamination. Also, extracts were collected into polypropylene conical tubes shown to be PFAS-free previously.

## Extractions and Soak Tests

To confirm that the EDGE was PFAS-free, a series of 12 sequential extractions were done to clean the EDGE. In these extractions, no Q-Cups were used and the EDGE was run according to the EDGE Method for EDGE System Analysis below. The 12th and final extract was used for analysis of the background contamination of PFAS on the EDGE by Pace Analytical.

## EDGE Method

Q-Disc: None\*

### Cycle 1

Extraction Solvent: Methanol  
Top Add: 15 mL  
Bottom Add: 5 mL  
Rinse: 0 mL  
Temperature: 65 °C  
Hold Time: 02:00 (mm:ss)

### Wash

Wash Solvent: Methanol  
Wash Volume: 10 mL  
Temperature: 65 °C  
Hold: 00:15 (mm:ss)

Individual samples of PEEK tubing, polypropylene tubing, Q-Discs, Q-Matrix Hydra, and Q-Cups were screened for potential contamination via soak tests. For five days, 0.5 g each of PEEK and polypropylene tubing, six G1 Q-Discs, six C9 Q-Discs, and six Q-Disc packaging separators were soaked in HPLC-grade methanol that was previously shown to be PFAS-free. A 5 g sample of Q-Matrix Hydra was soaked similarly. The soaking methanol for each sample type was subjected to analysis by Pace Analytical.

The Q-Cups were extracted in the EDGE without a Q-Disc using a two cycle method with two washes. The EDGE Method for Q-Cup Analysis parameters are shown below. The extract was subjected to analysis for its PFAS content by Pace Analytical.

\*We chose not to run a Q-Disc for this study however, in the software, it is required that you choose a Q-Disc.

## EDGE Method

Q-Disc: None\*

### Cycle 1

Extraction Solvent: Methanol/Water (80:20) with 0.3% ammonium hydroxide  
Top Add: 10 mL  
Bottom Add: 0 mL  
Rinse: 0 mL  
Temperature: 65 °C  
Hold Time: 03:00 (mm:ss)

### Cycle 2

Extraction Solvent: Methanol/Water (80:20) with 0.3% ammonium hydroxide  
Top Add: 10 mL  
Bottom Add: 0 mL  
Rinse: 0 mL  
Temperature: 65 °C  
Hold Time: 04:00 (mm:ss)

### Wash 1

Wash Solvent: Methanol  
Wash Volume: 10 mL  
Temperature: 50 °C  
Hold: 00:03 (mm:ss)

### Wash 2

Wash Solvent: Methanol/Water (80:20) with 0.3% ammonium hydroxide  
Wash Volume: 10 mL  
Temperature: - - -  
Hold: - :- -

## Results and Discussion

The 12th and final extract obtained from the EDGE system was analyzed. The list of compounds assessed, their detection limits, and the PFAS content from the extraction of the system and individual soak tests are in **Table 1** (page3). The data showed that the system was free of PFAS molecules below the level of detection after modification. Furthermore, soak tests of the samples of PEEK and polypropylene tubing, the G1 Q-Disc and C9 Q-Disc alone and together, their packaging separators, and Q-Matrix Hydra indicated that these consumables did not have any background PFAS contamination. Q-Cups were also shown to have no PFAS contamination. Thus, the EDGE is PFAS contamination-free.

## Conclusion

PFAS testing is extremely important because of the widespread nature of the compounds and their stability. Their use in the manufacture of materials that are commonly used in sample preparation and LC-MS analysis has presented a hurdle for their low level detection. Thus, instrumentation used in the extraction and analysis of PFAS molecules must be free of any prior PFAS contamination. The modifications made to the EDGE were effectively able to eliminate PFAS contamination. The consumables and Q-Cups used by the EDGE were also found to be PFAS free. The EDGE is an ideal choice for any laboratory seeking to automate their PFAS extractions without risking PFAS contamination.

**Table 1.** PFAS Compounds Assessed, Their Respective Detection Limit, and Soak Test Results

Compound	Limit of Detection (ng/L)	EDGE System	PEEK Tubing	Polypropylene Tubing	C9 Q-Disc	G1 Q-Disc	Q-Disc Separator	Q-Matrix Hydra	Q-Cup
9-chlorohexadecafluoro-3-oxanone-1-sulfonic acid (9Cl-PF3ONS)	25	ND	ND	ND	ND	ND	ND	ND	ND
11-chloroicosadecafluoro-3-oxaundecane-1-sulfonic acid (11Cl-PF30UdS)	25	ND	ND	ND	ND	ND	ND	ND	ND
1H, 1H, 2H, 2H-perfluorodecane sulfonic acid (8:2 FTS)	50	ND	ND	ND	ND	ND	ND	ND	ND
1H, 1H, 2H, 2H-perfluorooctane sulfonic acid (6:2 FTS)	50	ND	ND	ND	ND	ND	ND	ND	ND
1H,1H,2H,2H-perfluorododecane sulfonic acid (10:2 FTS)	50	ND	ND	ND	ND	ND	ND	ND	ND
1H,1H,2H,2H-perfluorohexane sulfonic acid (4:2 FTS)	50	ND	ND	ND	ND	ND	ND	ND	ND
Hexafluoropropylene oxide dimer acid (GenX)	50	ND	ND	ND	ND	ND	ND	ND	ND
4,8-dioxa-3H-perfluorononanoic acid (ADONA)	50	ND	ND	ND	ND	ND	ND	ND	ND
N-ethylperfluoro-1-octanesulfonamide (EtFOSA)	50	ND	ND	ND	ND	ND	ND	ND	ND
N-ethylperfluoro-1-octanesulfonamidoacetic acid (EtFOSAA)	50	ND	ND	ND	ND	ND	ND	ND	ND
2-N-ethylperfluoro-1-octanesulfonamido-ethanol (EtFOSE)	50	ND	ND	ND	ND	ND	ND	ND	ND
N-methylperfluoro-1-octanesulfonamide (MeFOSA)	100	ND	ND	ND	ND	ND	ND	ND	ND
N-methylperfluoro-1-octanesulfonamidoacetic acid (MeFOSAA)	50	ND	ND	ND	ND	ND	ND	ND	ND
2-N-methylperfluoro-1-octanesulfonamido-ethanol (MeFOSE)	50	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-1-butanefulfonic acid (PFBS)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-1-decanesulfonic acid (PFDS)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-1-heptanesulfonic acid (PFHpS)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-1-nonanesulfonic acid (PFNS)	53	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-1-octanesulfonamide (PFOSA)	50	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-1-pentanesulfonic acid (PFPeS)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluorohexanesulfonic acid (PFHxS)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-n-butanoic acid (PFBA)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-n-decanoic acid (PFDA)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-n-dodecanoic acid (PFDoA)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-n-heptanoic acid (PFHpA)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-n-hexanoic acid (PFHxA)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-n-nonanoic acid (PFNA)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-n-octanoic acid (PFOA)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-n-pentanoic acid (PFPeA)	25	ND	ND	ND	ND	ND	ND	ND	ND



Compound	Limit of Detection (ng/L)	EDGE System	PEEK Tubing	Polypropylene Tubing	C9 Q-Disc	G1 Q-Disc	Q-Disc Separator	Q-Matrix Hydra	Q-Cup
Perfluoro-n-tetradecanoic acid (PFTeDA)	50	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-n-tridecanoic acid (PFTrDA)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-n-undecanoic acid (PFUdA)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluorooctanesulfonic acid (PFOS)	38	ND	ND	ND	ND	ND	ND	ND	ND

**United States  
(Headquarters)**

800-726-3331  
704-821-7015  
info@cem.com

**France**

33 (01) 69 35 57 80  
info.fr@cem.com

**Germany, Austria,  
Switzerland**

(49) 2842-9644-0  
info@cem.de

**Ireland**

+353 (0) 1 885 1752  
info.ireland@cem.com

**Italy**

(39) 35-896224  
info.srl@cem.com

**Japan**

+81-3-5793-8542  
info@cemjapan.co.jp

**United Kingdom**

(44) 1280-822873  
info.uk@cem.com

**[www.cem.com](http://www.cem.com)**

© 2020 CEM Corporation  
All rights reserved. This may not be reproduced or published without written permission from CEM.