

# Extraction of Semi-Volatile Organic Compounds from Soil



## Abstract

Soil is one of the most common matrices in which semi-volatile organic compounds (SVOCs) are present. The extraction of these compounds from soil can be a lengthy and tedious process. The EDGE® is a revolutionary simple system for the rapid extraction of SVOCs from soil that is more than six times faster than other automated techniques. The Q-Cup® technology combines the process of Pressurized Fluid Extraction (PFE) and Dispersive Solid Phase Extraction (dSPE) in a single instrument that yields a fast, simple, and efficient extract. While maintaining all properties of PFE and adhering to US EPA 3545A, the EDGE also offers the benefits of the easy-to-assemble Q-Cup sample holder. The Q-Cup, with its unique open cell concept, creates a dispersive effect and promotes rapid extraction and filtration.

## Introduction

Semi-volatile organic compounds (SVOCs) are a subgroup of Volatile Organic Compounds (VOCs) that have high molecular weights and high boiling points. They are primarily composed of pesticides and herbicides. Prolonged exposure to these compounds, especially indoors, raises a public health concern, including several listed by the US EPA as hazardous air pollutants (HAPs). This classification applies to pollutants that can cause serious health effects, such as allergies, asthma, endocrine and thyroid disruption, reproductive toxicity, fetal and child development delays, and even cancer.

SVOCs have a broad set of chemical properties and structural features. These differences make it challenging to efficiently extract all analytes of interest with one method. Furthermore, the soil matrix from which the SVOCs are to be extracted often include multiple components, adding to the complexity of extraction. The EDGE with Q-Cup technology can effectively extract a difficult set of analytes from complex matrices with one simple method.

Traditional methods, such as Soxhlet, are very time consuming and use a large amount of solvent. Automated methods often require tedious sample preparation with complex sample holders. The EDGE is the fastest extraction system available and uses a minimal amount of solvent. The Q-Cup sample holder is comprised of two easy-to-assemble pieces, allowing the sample to be prepared in seconds.

Because of their persistent nature, SVOCs continue to accumulate and concentrate in our environment. To ensure our safety, these compounds need to be extracted, identified, and quantified. The accuracy of SVOC analysis is dependent on an efficient extraction. The EDGE yields an efficient extract that is filtered, cooled, and ready for analysis in under 10 minutes. US EPA 3545A is a method for the extraction of water-insoluble or slightly water-soluble volatile and semi-volatile compounds from soils, clays, sediments, sludges, and waste solids. The EDGE meets the requirements of US EPA 3545A and is preprogrammed for that method.

## Materials and Method

### Reagents

Sand, loam, and clay, purchased from Sigma Aldrich were spiked with SPEX CertiPrep TCLP Base/Neutral/Acid Extractable Spike Solution. The spike mix was also used to construct calibration curves. CRM 110-100 was purchased from Sigma Aldrich. Hexane/Acetone (1:1) was used for extraction rinse and wash.

### Sample Preparation

A portion of 15 grams of sand, loam, or clay spiked with 250 µl of spike solution was weighed into an assembled Q-Cup containing a S1 Q-Disc® stack (C9+G1+C9 sandwich). A 15 g portion of CRM 110-100 from Sigma Aldrich was also weighed into an assembled Q-Cup containing a S1 Q-Disc stack. The Q-Cups were placed in the EDGE removable rack, each with a collection vial, and the rack was positioned on the EDGE. The CEM method was started.

### EDGE Method for SVOCs from Soil

Q-Disc: S1 Q-Disc Stack (C9+G1+C9 sandwich)

#### Cycle 1

Extraction Solvent: Hexane/Acetone (1:1)  
Top Add: 20 mL  
Bottom Add: 10 mL  
Rinse: 0 mL  
Temperature: 100 °C  
Hold time: 03:00 (mm:ss)

#### Cycle 2

Extraction Solvent: Hexane/Acetone (1:1)  
Top Add: 20 mL  
Bottom Add: 10 mL  
Rinse: 0 mL  
Temperature: 100 °C  
Hold time: 03:00 (mm:ss)

#### Wash 1

Wash Solvent: Hexane/Acetone (1:1)  
Wash Volume: 30 mL  
Temperature: 120 °C  
Hold: 00:15 (mm:ss)

#### Wash 2

Wash Solvent: Hexane/Acetone (1:1)  
Wash Volume: 15 mL  
Temperature: - - -  
Hold: - :- -

### Analysis

The extract samples were concentrated to a volume of 10 mL. An aliquot of each extract was injected into the Agilent 7890A with a 5975C MSD for analysis, adhering to US EPA 8270. A Phenomenex ZB-5MSplus 30 m, 0.25 mm column was used.

## Results and Discussion

The EDGE extracted sand, loam, and clay samples in under 10 minutes, including filtration, cooling, and system washing. No post clean up or solvent exchange was necessary. The extracts were concentrated and injected into a GCMS for analysis.

**Figure 1** (page 3) is a representative GCMS chromatogram of a clay extraction showing clean separation of the SVOCs. The SVOCs analyzed are known to be difficult to extract, however the EDGE was able to achieve good extraction of all compounds in different soil types with a single method. The recovery data were measured via a 6-point calibration curve for each SVOC. The absolute recoveries for the SVOCs spiked in all three types of soil, sand, loam, and clay were acceptable. **Table 1** (page 3) shows the percent recoveries from spiked soil of some difficult-to-extract semi-volatile organic compounds. The extraction of CRM 110-100 via EDGE was comparable to Soxhlet as seen in **Table 2** (page 3).

## Conclusion


To be in the acceptable recovery range, the percent recovery data for the extraction of SVOCs from soil should be within 80-120%. The EDGE utilized its Q-Cup technology to yield acceptable recoveries for the extraction of SVOCs from three different types of soil and a CRM. Furthermore, the EDGE method used less solvent and time than traditional extraction techniques such as Soxhlet. Lastly, since the EDGE is an automated system, it eliminates the human error often found with other extraction techniques. The EDGE was able to efficiently and accurately extract SVOCs from soil.

**Table 1.** Percent Recovery Data for 15 g Spiked Sand, Loam, and Clay

Compound	Sand	Clay	Soil
Pyridine	100	88	93
1,4-dichlorobenzene	88	88	96
2-methylphenol	84	95	115
3-methylphenol	90	102	104
Hexachlorobutadiene	86	92	97
2,4,6-trichlorophenol	90	105	103
2,4,5-trichlorophenol	89	113	99
2,4-dinitrotoluene	90	102	99
Hexachlorobenzene	86	86	81

**Table 2.** Percent Recovery Data as Compared to Soxhlet for 15 g CRM 110-100

Compound	% Soxhlet
2-nitroaniline	94
2,4-dinitrotoluene	114
Dibenzofuran	92
Fluorene	105
Bis-2-ethylhexyl-phthalate	97



**Figure 1.** GCMS Chromatogram of the SVOCs Extracted from Clay

**United States  
(Headquarters)**

800-726-3331  
704-821-7015  
Fax: 704-821-7894  
info@cem.com

**France**

33 (01) 69 35 57 80  
Fax: 33 (01) 60 19 64 91  
info.fr@cem.com

**Germany, Austria,  
Switzerland**

(49) 2842-9644-0  
Fax: (49) 2842-9644-11  
info@cem.de

**Ireland**

+353 (0) 1 885 1752  
Fax: +353 (0) 1 885 1601  
info.ireland@cem.com

**Italy**

(39) 35-896224  
Fax: (39) 35-891661  
info.srl@cem.com

**Japan**

+81-3-5793-8542  
Fax: +81-3-5793-8543  
info@cemjapan.co.jp

**United Kingdom**

(44) 1280-822873  
Fax: (44) 1280-822873  
info.uk@cem.com

**[www.cem.com](http://www.cem.com)**

© 2020 CEM Corporation  
All rights reserved. This may not be reproduced or published without written permission from CEM.