


Automated Microwave-Enhanced Synthesis of Glycopeptides with O-Linked Glycans



Summary

Glycopeptides containing O-linked glycans can be synthesized quickly with good purity using microwave-enhanced SPPS. Preparation of a peptide substrate of Wbwa1 (a sialic acid glycosyltransferase), was completed in 3 h 10 min with 68% purity on the Liberty Blue™. Synthesis of an analog of APF (antiproliferative factor)² conjugated to the TUS nuclear localization sequence³ was achieved in 1 h 8 min with 73% purity using the Liberty PRIME™.

Introduction

Glycosylation is an important post-translational and co-translational modification that involves the covalent modification of proteins or peptides with mono, di, or oligosaccharides.⁴ Glycoproteins have diverse biological functions ranging from structural proteins and hormones (like collagen and thyroid-stimulating hormone) to enzymes and proteins involved in regulating immune response (such as various phosphatases and antibodies).⁵ As the critical role of glycosylation becomes further illuminated, interest in the chemical synthesis of glycopeptides continues to grow steadily. Peptides containing O-glycosylated amino acids can be prepared by SPPS by using Fmoc protected Ser or Thr (**Figure 1**). Application of microwave energy to the synthesis of glycopeptides containing O-linked glycans allows for more efficient coupling which leads to rapid synthesis times and high purity (CarboMAX™).⁶




Figure 1: Fmoc-(S)-Ser(GlcNAc(Ac)₃-β-D)-OH (left) and Fmoc-(S)-Thr(GlcNAc(Ac)₃-β-D)-OH (right)

Materials and Methods

Reagents

The following Fmoc amino acids were obtained from CEM Corporation (Matthews, NC) and contain the indicated side chain protecting groups: Ala, Arg(Pbf), Asp(OMpe), Gly, His(Boc), Ile, Leu, Lys(Boc), Pro, Ser(tBu), Thr(tBu), and Val. Rink Amide ProTide™ LL resin and Fmoc-Ala-Wang LL resin were also obtained from CEM Corporation. Fmoc-Ser(GlcNAc(Ac)₃-β-D)-OH and Fmoc-Thr(GlcNAc(Ac)₃-β-D)-OH were purchased from EMD Millipore (Burlington, MA). N,N'-Diisopropylcarbodiimide (DIC), piperidine, pyrrolidine, trifluoroacetic acid (TFA), 3,6-dioxo-1,8-octanedithiol (DOTD), and triisopropylsilane (TIS) were obtained from Sigma-Aldrich (St. Louis, MO). Dichloromethane (DCM), N,N-Dimethylformamide (DMF), anhydrous diethyl ether (Et₂O), and acetic acid were obtained from VWR (West Chester, PA). LC-MS grade water (H₂O) and LC-MS grade acetonitrile (MeCN) were obtained from Fisher Scientific (Waltham, MA).

Peptide Synthesis: WbwA Substrate, PAHGV-Thr(GlcNAc(Ac)₃-β-D)-SAPDTRPAPGSTAP-NH₂

The peptide was synthesized on a 0.05 mmol scale using the CEM Liberty Blue automated microwave peptide synthesizer on Rink Amide ProTide resin LL (0.18 meq/g substitution). Deprotection was performed with piperidine in DMF. Coupling reactions were performed in 5-fold excess of Fmoc-AA with DIC and Oxyma Pure (with DIEA) in DMF (CarboMAX).⁶ Fmoc-Thr(GlcNAc(Ac)₃-β-D)-OH was used for the glycosylated amino acid. Cleavage was performed at room temperature for 2 h using TFA/H₂O/TIS/DODT. Following cleavage, the peptide was precipitated with Et₂O and lyophilized overnight.

Peptide Synthesis: TUS-APF Analog, KLKIKRPVK-Ser(GlcNAc(Ac)₃-β-D)-VPAVVVA-CO₂H


The peptide was synthesized on a 0.10 mmol scale using the CEM Liberty PRIME automated microwave peptide synthesizer on Fmoc-Ala-Wang LL resin (0.31 meq/g substitution). Deprotection was performed with pyrrolidine in DMF. Coupling reactions were performed in 5-fold excess of Fmoc-AA with DIC and Oxyma Pure in DMF (CarboMAX).⁶ Fmoc-Ser(GlcNAc(Ac)₃-β-D)-OH was used for the glycosylated amino acid. Cleavage was performed at room temperature for 2 h using TFA/H₂O/TIS/DODT. Following cleavage, the peptide was precipitated with Et₂O and lyophilized overnight.

Peptide Analysis


The peptides were analyzed on a Waters Acquity UPLC system with PDA detector equipped with an Acquity UPLC BEH C8 column (1.7 mm and 2.1 x 100 mm). The UPLC system was connected to a Waters 3100 Single Quad MS for structural determination. Peak analysis was achieved on Waters MassLynx software. Separations were performed with a gradient elution of 0.05% TFA in (i) H₂O and (ii) MeCN.

Results

Microwave-enhanced SPPS of WbwA substrate on the Liberty Blue automated microwave peptide synthesizer produced the target peptide in 68% purity (**Figure 2**).

**Figure 2:** UPLC Chromatogram of WbwA Substrate

Microwave-enhanced SPPS of TUS-APF analog on the Liberty PRIME automated microwave peptide synthesizer produced the target peptide in 73% purity (**Figure 3**).

**Figure 3:** UPLC Chromatogram of TUS-APF analog**Conclusion**

Glycopeptides with O-linked glycans can be synthesized rapidly and efficiently using automated microwave-enhanced SPPS. A WbwA peptide substrate was synthesized in 3 h 10 min with 68% purity using the Liberty Blue peptide synthesizer. On the Liberty PRIME, microwave-enhanced SPPS affords a TUS-APF analog peptide in 1 h 8 min with a purity of 73%. Additionally, this method reduces the overall synthesis cost for O-linked glycopeptides by using less than 2 equivalents of the very expensive Fmoc-protected glycoamino acids.

References

- (1) Czuchry, D.; Desormeaux, P.; Stuart, M.; Jarvis, D. L.; Matta, K. L.; Szarek, W. A.; Brockhausen, I. *J. Bacteriol.* **2015**, *197* (24), 3760–3768.
- (2) Siddiqui, M. A.; Ambre, S.; Keay, S. K.; Rhyne, J. M.; Zhang, C.-O.; Barchi, J. *J. ACS Omega* **2017**, *2* (9), 5618–5632.

(3) Kaczmarczyk, S. J.; Sitaraman, K.; Hill, T.; Hartley, J. L.; Chatterjee, D. K. *PLoS One* **2010**, 5 (1), e8889.

(4) Spiro, R. G. *Glycobiology* **2002**, 12 (4), 43R–56R.

(5) Murray, R. K.; Granner, D. K.; Rodwell, V. W. *Harper's Illustrated Biochemistry*, 27th ed.; McGraw-Hill Medical, 2006.

(6) CEM Application Note (APO124) - "CarboMAX - Enhanced Peptide Coupling at Elevated Temperature."

**United States
(Headquarters)**

800-726-3331
704-821-7015
Fax: 704-821-7894
info@cem.com

France

33 (01) 69 35 57 80
Fax: 33 (01) 60 19 64 91
info.fr@cem.com

**Germany, Austria,
Switzerland**

(49) 2842-9644-0
Fax: (49) 2842-9644-11
info@cem.de

Ireland

+353 (0) 1 885 1752
Fax: +353 (0) 1 885 1601
info.ireland@cem.com

Italy

(39) 35-896224
Fax: (39) 35-891661
info.srl@cem.com

Japan

+81-3-5793-8542
Fax: +81-3-5793-8543
info@cemjapan.co.jp

United Kingdom

(44) 1280-822873
Fax: (44) 1280-822873
info.uk@cem.com

www.cem.com

© 2019 CEM Corporation
All rights reserved. This may not be
reproduced or published without
written permission from CEM.