



Step 1 of 2

**Procedure**

Weigh 0.1 g of the sample into the digestion vessel. Add 3 mL of HNO<sub>3</sub> and 5 mL of HF. Gently swirl the mixture and wait approximately 15 minutes before closing the vessel.

**Notes**

This procedure uses hydrofluoric acid. If it is necessary to complex the residual hydrofluoric acid or redissolve insoluble fluorides formed, an additional complexation step with boric acid should be used. This procedure can be found in the One Touch Method note entitled “Boric Acid HF Neutralization”.

**Recommended Equipment**

MARS 6 iWave

**Recommended Vessels**

iPrep

**Reagents**

HNO<sub>3</sub>  
HF

**Max Sample Weight**

0.1 g

**Sample Type**

Organic

**Control Type**

Ramp to Temperature

**Method Type**

One Touch

**Heating Program**

Stage	Temp (°C)	*Ramp (mm:ss)	Hold (mm:ss)	Pressure (psi)	* Power (W)	Stirring
1	210	25:00	25:00	N/A	700-1800	Off

\* Ramp times and power may vary depending on the type and number of vessels.



Step 2 of 2

**Procedure**

Allow vessel to cool. Add 2.5 g H<sub>3</sub>BO<sub>3</sub> and 25 mL of deionized H<sub>2</sub>O into the vessel that contains the sample and acid.

**Notes**

This procedure can be used if it is necessary to complex the residual hydrofluoric acid or redissolve insoluble fluorides formed by reaction of certain analytes with hydrofluoric acid.

**Reagents**

H<sub>3</sub>BO<sub>3</sub>  
DI H<sub>2</sub>O

**Sample Type**

Organic

**Control Type**

Ramp to Temperature

**Method Type**

One Touch

**Heating Program**

Stage	Temp (°C)	*Ramp (mm:ss)	Hold (mm:ss)	Pressure (psi)	* Power (W)	Stirring
1	170	25:00	15:00	N/A	700-1800	Off

\* Ramp times and power may vary depending on the type and number of vessels.

**Results**

Sample was clear and particle free upon dilution to 50 mL.

**General Precaution**

- a) If using HF, follow restrictions listed in HF Addendum.
- b) This procedure is a reference point for sample digestion using a CEM system and may need to be modified or changed to obtain the required results on your sample.
- c) If using a vessel other than the recommended choice, adjust sample size and pressure limit to values appropriate for the vessel chosen.
- d) Manual venting of CEM vessels should be performed when wearing hand/eye/body protection and when the vessel contents are at or below room temperature to avoid the potential for chemical burns. Always point the vent hole away from the operator.
- e) If programming as one touch, the ramp time and power will be automatically determined based on the number and type of vessels detected.