

Part Numbers: S001-C, S001-D

1. IDENTIFICATION OF PRODUCT AND THE COMPANY

1.1 Product identifiers

Product name: Oxyma Pure; Ethyl (hydroxyimino)cyanoacetate Chemical Name: Ethyl (hydroxyimino)cyanoacetate

1.2 Details of the supplier of the safety data sheet

Company: CEM Corporation 3100 Smith Farm Rd, Matthews, NC 28104

Telephone: 704-821-7015 Fax: 704.821.7894 24 hour access: +1-704-564-4503 E-mail address: peptide.support@cem.com

2. HAZARDS IDENTIFICATION

- **2.1 Classification of the substance or mixture** Not a hazardous substance or mixture.
- **2.2 GHS Label elements, including precautionary statements** Not a hazardous substance or mixture.
- 2.3 Hazards not otherwise classified (HNOC) or not covered by GHS None

3. COMPOSITION/INFORMATION ON INGREDIENTS

3.1 Substances:

Synonyms: Oxyma Pure Formula: $C_5H_6N_2O_3$ Molecular Weight: 142.11 g/mol CAS No.: 3849-21-6 EC-No. : 223-351-3

No ingredients are hazardous according to OSHA criteria. No components need to be disclosed according to the applicable regulations.

4. FIRST AID MEASURES

4.1 Description of first aid measures

If inhaled:

If breathed in, move person to fresh air. If not breathing, give artificial respiration.

In case of skin contact Wash off with soap and plenty of water. In case of eye contact Flush eyes with water as a precaution. If swallowed Never give anything by mouth to an unconscious person. Rinse mouth with water.

4.2 Most important symptoms and effects, both acute and delayed

The most important known symptoms and effects are described in the labelling (see section 2.2) and/or in section 11

4.3 Indication of immediate medical attention and special treatment needed No data available.

5. FIRE-FIGHTING MEASURES

5.1 Extinguishing media

Suitable extinguishing media Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide

5.2 Special hazards arising from the substance or mixture Carbon oxides, nitrogen oxides (NOx)

5.3 Advice for firefighters

Wear self-contained breathing apparatus for firefighting if necessary.

5.4 Further information

No data available

6. ACCIDENTAL RELEASE MEASURES

6.1 Personal precautions, protective equipment, and emergency procedures Avoid dust formation. Avoid breathing vapors, mist or gas.

6.2 Environment precautions

Do not let product enter drains.

6.3 Methods and materials for containment and cleaning up Sweep up and shovel. Keep in suitable, closed containers for disposal.

6.4 Reference to other sections

For disposal see section 13.

7. HANDLING AND STORAGE

7.1 Precautions for safe handling

Provide appropriate exhaust ventilation at places where dust is formed. Normal measures for preventive fire protection. For precautions see section 2.2

7.2 Conditions for safe storage, including any incompatibilities

Keep container tightly closed in a dry and well-ventilated place.

7.3 Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific uses are stipulated

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

8.1 Control parameters

Components with workplace control parameters

Contains no substances with occupational exposure limit values.

8.2 Exposure controls

Appropriate engineering controls General industrial hygiene practice.

Personal protection equipment

Eye/face protection

Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

Skin protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands.

Body protection

Choose body protection in relation to its type, to the concentration and amount of dangerous substances, and to the specific work-place. The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

Respiratory protection

Respiratory protection is not required. Where protection from nuisance levels of dusts are desired, use type N95 (US) or type P1 (EN 143) dust masks. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Control of environmental exposure

Do not let product enter drains.

9. PHYSICAL AND CHEMICAL PROPERTIES

9.1 Information on basic physical and chemical properties

a) Appearance	Form: crystalline	
	Color: beige	
b) Odor	No data available	
c) Odour Threshold	No data available	
d) pH value	No data available	
e) Melting point/freezing point		
	Melting point/range: 130 - 132 °C (266 - 270 °F)	
f) Initial boiling point and b	oiling range	
	No data available	

g) Flash point h) Evapouration rate	No data available No data available	
i) Flammability	No data available	
j) Upper/lower flammability or explosive lim		
	No data available	
k) Vapor pressure	No data available	
l) Vapour density	No data available	
m) Relative density	No data available	
n) Water solubility	No data available	
 o) Partition coefficient n-octanol/Water 		
	No data available	
p) Auto-ignition Temperature		
	No data available	
q) Decomposition Temperature		
	No data available	
r) Viscosity	No data available	
s) Explosive properties	No data available	
t) Oxidizing properties	No data available	

9.2 Other safety information

No data available

10. STABILITY AND REACTIVITY

10.1 Reactivity

No data available

10.2 Chemical stability

Stable under recommended storage conditions

10.3 Possibility of hazardous reactions

No data available

- **10.4 Conditions to avoid** No data available
- 10.5 Material to avoid Strong acids, Strong bases

10.6 Hazardous decomposition procedures:

Other decomposition products - no data available In the event of fire: see section 5

11. TOXICOLOGICAL INFORMATION

11.1 Information on toxicological effects Acute toxicity No data available

Inhalation: No data available

Dermal: No data available

Skin corrosion/irritation No data available

Serious eye damage/eye irritation No data available

Respiratory or skin sensitisation No data available

Germ cell mutagenicity No data available

Carcinogenicity

- IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.
- ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.
- NTP: No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.
- OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.

Reproductive toxicity

No data available

Specific target organ toxicity - single exposure No data available

Specific target organ toxicity - repeated exposure No data available

Aspiration hazard No data available

Additional Information RTECS: Not available

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

12. ECOLOGICAL INFORMATION

- 12.1 Toxicity No data available
- **12.2 Persistence and degradability** No data available
- 12.3 Bioaccumulative potential No data available

12.4 Mobility in soil

No data available

12.5 Results of PBT and vPvB assessment PBT/vPvB assessment not available as chemical safety assessment not required/not conducted

12.6 Other adverse effects

No data available

13. DISPOSAL CONSIDERATIONS

13.1 Waste treatment methods

Product

Offer surplus and non-recyclable solutions to a licensed disposal company.

Contaminated packaging

Dispose of as unused product.

14. TRANSPORT INFORMATION

DOT (US) Not dangerous goods

IMDG Not dangerous goods

IATA Not dangerous goods

15. REGULATORY INFORMATION

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

No SARA Hazards

16. OTHER INFORMATION

Further information

Copyright 2015 CEM Corporation, is granted to make unlimited paper copies for internal use only. The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. The information in this document is based on the present state of our knowledge and is applicable to the product with regard to appropriate safety precautions. It does not represent any guarantee of the properties of the product. CEM Corporation shall not be held liable for any damage resulting from handling or from contact with the above product.