



Systems

ORACLE

Equipment

Analytical Balance
CEM Square Glass Fiber Pads
CEM Trac Film
Heater Block
ORACLE Station

ORACLE Parameters

QuickPrep™ Temp	QuickPrep™ Time	Warm Delay
70	30	5

Sample Preparation

Procedure

On external balance, tare glass fiber sample pad then add 1.9-2.1 grams of sample to the pad. Prepare the samples using the ORACLE Station, input the weight of sample on the "run screen" and drop sample into QuikPrep™.