

## Microwave Extraction System







## Better PFAS Extraction at Your Fingertips

High Throughput • User-friendly Software • Scientist Centered Design

The MARS 6 PFAS Extraction system is specifically designed with PFAS extraction in mind. Process up to 40 samples simultaneously, yielding the highest throughput of any system available. The simple turntable utilizes 50 or 15 mL PFAS-free polypropylene centrifuge tubes. Effortless operation facilitates transfer from sample to SPE or directly to analysis. PFAS determination has never been easier.

## Simple Process

- 1. Select a Method.
- 2. Matrix is Extracted.
- 3. Sample is Cooled.

- High Throughput System 40 samples in a single batch
- Disposable vessels for easy cleaning
- Meets the requirements of draft method EPA 1633
- All components in contact with sample are tested to be PFAS-free
- Compatible with a wide variety of sample types



Visit cem.com/mars-6-pfas for more details