





	Liberty Lite™	Liberty Blue™	Liberty Blue HT12™	Liberty PRIME™
				
<b>Cycle Time</b>	15 minutes	4 minutes	4 minutes	2 minutes
<b>Waste/Cycle</b>	40 mL (0.1 mmol)	16 mL (0.1 mmol)	16 mL (0.1 mmol)	8.5 mL (0.1 mmol)
<b>Scale Range</b>	0.005 - 5 mmol	0.005 - 5 mmol	0.005 - 5 mmol	0.005 - 5 mmol
<b>Peptide Positions</b>	1	1	12	24
<b>Amino Acid Positions</b>	20	27	27	27
<b>Other Positions</b>	4	4	4	4
<b>Fluid Delivery</b>	<b>Flex-Add™:</b> Amino Acids & Activators  <b>Timed Delivery:</b> Wash & Deprotection	<b>Flex-Add:</b> Amino Acids & Activators  <b>Timed Delivery:</b> Wash & Deprotection	<b>Flex-Add:</b> Amino Acids & Activators  <b>Timed Delivery:</b> Wash & Deprotection	<b>Flex-Add:</b> Amino Acids  <b>Pre-Calibrated Pumps:</b> Wash, Deprotection, & Oxyma
<b>Dimensions</b>	· 20" W x 18" D x 30" H · 51 cm x 46 cm x 76 cm	· 20" W x 18" D x 30" H · 51 cm x 46 cm x 76 cm	· 27" W x 18" D x 30" H · 69 cm x 46 cm x 76 cm	· 42" W x 18" D x 30" H · 107 cm x 46 cm x 76 cm
<b>Upgrades</b>	Upgrade to the Liberty Blue	· HT12 · HT24	· HT24	N/A
<b>Accessories</b>	Razor Cleavage System	· Integrated Camera, · Razor Cleavage System · Flex-Add Large Scale	· Integrated Camera · Razor Cleavage System · Flex-Add Large Scale	· Integrated Camera · Razor Cleavage System



### Razor®

## Cleave peptides in only 30 minutes.

The Razor is the perfect compliment to a Liberty Peptide Synthesizer. Cleave peptides in only 30 minutes at elevated temperature while protecting your synthesizer from volatile TFA fumes.